BCI practical course: "Hello World" & ERP Viewer

Jason Farquhar

<J.Farquhar@donders.ru.nl>

Learning Goals

Understand:

- What is needed to make a BCI, i.e. progress tracking, data acquisition, annotation and processing, stimulus presentation, and an overall process scheduler/sequencer
- How to used event-driven programming ideas coupled to a global shared event pool (blackboard) to provide these facilities
- How the fieldtrip buffer provides the event blackboard which is used for inter-process communication.

Know how to:

- What the struct of an 'event' is and how to use it to annotate data with experiment relevant event information
- present simple visual stimulus/feedback to the user/experimenter
- How to wait for specific events, get the necessary data, process it and post the updated results back to the event blackboard
- Startup the buffer and an experiment control Matlab, and how to connect these processes to provide a basic BCI
- Test your experiment with simulated data generated by the signal-proxy
- Debug your experiment when it fails!



Today's Plan

- Discussion: What do we need to make a BCI?
- Introduction to the Buffer-BCI framework

break

- Hands-on 1: Hello World
- Hands-on 2: Sequenced Sentences

break

- Hands-on 3: Echo-server (event IPC)
- Hands-on 4: Visual ERP Viewer




Discussion: What do we need to make a BCI?

 Based on your prior knowledge and experience with the hands on demo we've did last time.


- Discuss: What do we require to make a BCI system?
- Think about:
 - Hardware requirements?
 - Software requirements?
 - Information flows?



BCI information flow




Gross structure of a typical BCI experiment



BCI terminology (our group!)

- Epoch/Trial
 - Single BCI prediction
 - e.g. 1-imagined movement, 1 visual-speller flash
- Sequence
 - Short group of epochs (~1min)
 - v. short breaks 1-2sec between epochs (usually automatic)
 - short (usually self-paced) subject break between sequences (~10sec)
- Block/Run/Phase
 - Short group of sequences (>10min)
 - long (~1-2min) subject break between blocks
 - e.g. cap-fitting, calibration, classifier training, on-line use
- Session during one cap-fitting
- Experiment imagined movement, visual-speller etc.

Flow chart of an individual epoch in a simple BCI experiment



Requirements: what do you need to build a BCI?

- 1) Way of **tracking** where we are in execution of the experiment flowchart, i.e. block, sequence, epoch number.
- 2)Way of **annotating** data to what the subject was experiencing/doing at that time with what was measured from their brain/body, e.g. LH movement, reading instruction, watching queue, etc.
- **3)**Data acquisition: Drivers to extract data from hardware (and combine data fro different hardware sources)
- **4)Stimulus Generation**: makes stimuli that the subject will experience, for subject instruction, feedback, event-related stimuli
- 5)Something to **process** the signals, firstly to train the classifier, and secondly to decode the users mental state, i.e. do the BCl bit ;-)
- 6) Scheduler (sequencer?) to tie it all these bits together,
 - so the correct functions, i.e. stimulus display, signal processing, are executed
 - at the correct position in the experiment flowchart
 - · based on the right bits of measured datacc



Summary

- To build a BCI we need a system to; track our progress through the experiment, acquire, annotate and process data, present the stimuli and schedule all these processes in an appropriate way.
- Next we introduce a Matlab based system which provides these facilities.

Buffer-BCI Framework


- •We can break the requirements into 4 largely independent communicating processes:
 - 1 Data-acquisitation & annotation
 - ·Get data from hardware
 - ·Attach annotations (markers, events) to particular data sample
 - 2)Experiment control (scheduling)
 - · Control the flow of the experiment
 - 3)Stimulus generation
 - · Make stimuli when requested by the expt controller
 - · Make feedback based on predictions generated by the sig-processor
 - 4)Signal Processing
 - Process the data based on the annotationss, and generate predictions



Buffer-BCI Framework

Basic Idea:

- set of independent processes
- any process can send/recieve dataannotation events
- events are visible to all other processes
- Processes communication implemented by sending recieving events
- (N.B. As all events are saved with the data, annotations are automatically archived for later off-line use).
- Similar in concept to that used in 'Blackboard architectures' for AI, see < en.wikipedia.org/wiki/Blackboard system>



Ft-buffer based Implementation

- Buffer-BCI framework implemented using the fieldtrip-buffer system (fieldtrip.fcdonders.nl/development/realtime)
- Ft-buffer provides:
 - Drivers for data-acquisation
 - 1)buffer storage for data (~last 1 minute data)
 - 2)buffer storage for events (~last 50 events)

Idea:

- Buffer events store represents the blackboard used for inter-process communication (IPC)
- Every event has timestamp (sample number) used for data-annotation




Buffer-BCI Framework



Basic Idea:

- Client/Server architecture
- set of independent processes
- any process can send/recieve dataannotation events
- events are visible to all other processes
- Processes communication implemented by sending receiving events
- (N.B. As all events are saved with the data, annotations are automatically archived for later off-line use).
- Similar in concept to that used in 'Blackboard architectures' for Al, see < en.wikipedia.org/wiki/Blackboard_system>







0

Ft-buffer based Implementation

- Buffer-BCI framework implemented using the fieldtrip-buffer system (fieldtrip.fcdonders.nl/development/realtime)
- Ft-buffer provides:
 - Drivers for data-acquisation
 - 1) buffer storage for data (~last 1 minute data)
 - 2) buffer storage for events (~last 50 events)


Idea:

- Buffer events store represents the blackboard used for interprocess communication (IPC)
- Every event has timestamp (sample number) used for dataannotation





Typical Structure of BCI system








client/server architecture

- 1)Buffer server: for data storage, annotation (IPC) and archive
- 2)Acquisition driver: to access Decoding (or simulate) the measurement (Sig-proc) hardware
- 3)BCI Application.

Normally, split into 2/3 pieces:

- 1) Signal-processing
- 2) Stimulus presentation
- 3) Experiment control (commonly part of stimulus-presentation)









Buffer_bci Quickstart




client/server architecture

- 1) Buffer server: for data storage, annotation (IPC) and archive
- 2) Acquisition driver: to access (or simulate) the measurement hardware
- 3) BCI Application.

Normally, split into 2/3 pieces:

- 1) Signal-processing
- 2) Stimulus presentation
- 3) Experiment control (commonly part of stimulus-presentation)







Running buffer-bci code:

You need to have (at least) the following processes running:

- 1) Buffer server
- 2) Data-acquisation system (real or simulated)
- 3) BCI Application

Quickstart: 1) Buffer Server

- Saving server (saves all data and events)
 dataAcq/startJavaBuffer.{sh/bat}
- Server (saves nothing): dataAcq/startJavaNoSaveBuffer.{sh/bat}

C-based implementations:

- Saving: dataAcq/startBuffer.{sh/bat}
- NoSaving: dataAcq/startNoSaveBuffer.{sh/bat}



Quickstart: 2) Data-acquisation

Simulated data: (for debugging/testing)

• Java-based: dataAcq/startJavaSignalProxy. {sh/bat} (Alternative implementions Matlab: startMatlabSignalProxy, c:startSignalProxy)

Real data: (for debugging/testing)

- FilePlayback: dataAcq/startJavaFileProxy. {sh/bat}
- Microphone: dataAcq/startAudio.{sh/bat}

EEG amplifier data: (for real usage)

- **Biosemi:** dataAcq/startBiosemi.{sh/bat}
- Mobita: dataAcq/startMobita. { sh/bat }
- OpenBCI: dataAcq/startOpenBCI.{sh/bat}
- Many others... (look in dataAcq/startXXXXX. {bat/sh})



Quickstart: 3) BCI-Application

- Generally what you will write yourself in this course;-)
- Typically in 2 (or more) separate processes:

1) Stimulus presentation:

- Controls what the user sees/hears, i.e. Ul
 - Show App, show feedback, control output-devices, etc.

2) Signal processing:

- Processes the EEG data and events to generate predictions for UI feedback.
 - Collect calibration data, train classifier, generate predictions

Quickstart Example: Games-demo

1) Buffer Server:

```
dataAcq/startJavaNoSaveBuffer.{sh,bat}
```

2) Data Acquisation (simulated):

```
dataAcq/startJavaSignalProxy.{sh,bat}
```

3) BCI Application(1) -- Stimulus Presentation:

```
games/runGame.{bat,sh}
```

4) BCI Application(2) – Signal Processing:

```
games/startSigProcBuffer.{bat,sh}
```

Useful (debugging) Functions:

 Seeing the what events are sent and when is important for debugging experiment flow.

[]=eventViewer(host,port,mtype,mval)

- Print all events matching (mtype, mval)
- (mtype,mval) are as used in matchEvents

N.B. you can run the eventViewer directly using: dataAcq/startEventViewer.{bat,sh}

EventViewer in different languages

Basic event viewer has been implemented in multiple languages

- MATLAB: utilities/eventViewer.m
- JAVA: java/echoClient/eventViewer.java
- Python: python/echoClient/eventViewer.py
- C:c/echoClient/eventViewer.c



Hands-on 1: Event Echo-Server

Events for IPC

- As well as being used for data annotation, Events are used for inter-processcommunication,
 - e.g. to communicate results from signal-processing to stimulus presentation
- To use events in this way, each process must
 - monitor for new events
 - filter out the events it should react to
 - send response events



(Key concept) event structure

sample

 time at which event occurred in samples from start of experiment

type

arbitrary event type (usually a string)

value

 arbitrary event value (usually string or number)

duration (optional)

duration of the event in samples

offset (optional)

- zero-time for the event.
- Usually, offset from sample at which the event actually started.

```
Examples:
Visual speller "flash";
ev=struct('sample',123,...
'type','stimulus.flash',...
'value',[0 0 1 0 0],...
'offset',0,'duration',0)
```

```
Classifier prediction:
ev=struct('sample',123,...
'type','prediction',...
'value',[-1 -1 -1 1 -1],...
'offset',0,'duration',0)
```

```
Imagined Movement event:
ev=struct('sample',123,...
'type','stimulus.move',...
'value','left-hand',...
'offset',0,'duration',300)
```

```
Compact notation: s:123,t:'stimulus.flash',v:[0 0 1 0 0],o:0,d:0
```

(key functions) Event manipulation

evt=mkEvent(type,value,[sample,offset,duration])

make a buffer event, with sensible defaults

sendEvent

- evt=sendEvent(type,value,[sample,offset,duration,host,port])
- evt=sendEvent(evt,[host,port])
- Send event to the buffer on machine host at port.

mi=matchEvents(*evts*,*mtype*,*mval*)

- Find events with type mtype and value mval in evts a vector of event structures.
- mtype can be cell-array of types to match, e.g. {'type1' 'type2'}
- mval can be cell-array of values to match, e.g. {'val1' 10 'val3'}
- match if any mtype matches and any mval matches,
 - i.e. above matches (t:'type1',v:10), (t:'type2',v:10),(t:'type1',v:'val1')
- mi is logical vector of which evts matched
- N.B. Empty ([]) or '*' mtype/mvalue matches everything



Hands-on 1: Event Echo-Server

Experiment Task

- Write a simple echo-server which:
 - Connects to the Buffer server
 - Waits for any incomming event, and
 - Responds by sending a 'echo' event with the same value but type='echo'
 - Quits if it recieves and event with type='exit'
- N.B. Don't 'echo' your own echo events!

Assignment:

• start from : echoServer skel.m

N.B. send 'keyboard' events by pressing keys in the signal proxy window.



Useful Functions:

- Setup MATLAB paths and connect to buffer:
 - See the header code block in echoServer_skel.m
 - This will initialize the paths:

```
run ../../utilities/initPaths.m
```

 Then try to connect to the buffer server until valid header is returned.

```
while (isempty(hdr) | ....
```

 Then initialize some utility functions for highprecision timing

```
initsleepSec();initgetwTime();
```

Useful Functions:

```
[devents,state]=...
buffer_newevents(host,port,state,mtype,mval,timeout_ms)
```

- wait for any new events matching (mtype, mval)
 - Matching done by matchEvents
 mtype can be cell-array of types to match, e.g. {'type1' 'type2'}
 mval can be cell-array of values to match, e.g. {'val1' 10 'val3'}
 match if any mtype matches and any mval matches
- return the matched events in the vector of structure(s) devents
- state is the match state, used to track which events have been processed between function calls
- Return after timeout_ms milliseconds even if no matching events found

Debug/Test your echoSever

- Start an event viewer: dataAcq/startJavaEventViewer. { sh, bat }
- Generate events to be echo'ed:
 - Directly from Matlab:
 - Start 2nd Matlab
 - Connect to buffer
 - Send events: sendEvent('type', 'value')
 - OR:

From the signal-proxy:

- Start the Matlab signal proxy:
 dataAcq/startMatlabSignalProxy.{sh,bat}

- Press keys in the SignalProxy window to generate {t:keyboard, v:key} events

Echo-Server in different languages

Basic echo-server example has been implemented in multiple languages

- MATLAB: echoClient/matlabclient.m
- JAVA: java/echoClient/javaclient.java
- C#:csharp/echoClient/csharpclient.cs
- Python: python/echoClient/pythonclient.py
- C:c/echoClient/cclient.c

Hands-on 2: "Hello World"

Experiment Task

- Display the string "Hello World" (or any other prespecified string) on the screen, and wait for a key to be pressed to exit
- Send events to annotate what has happened, e.g. startup, string display, key-pressed, shutdown etc.

Method:

- Start from the 'helloworld_skel.m' function skeleton
 - contains initialisation code to connect to the ft_buffer
 - Some examples of functions you may find useful



Note: event timestamps

- Accurate event time-stamps are critical for evoked potential analysis
 - >10ms event jitter causes significant reduction in signal quality
- However,
 - data-acquisation may only send data every >20ms
 - And this data may be subject to additional network delays of >20ms
- Stop this jitter reducing time-stamp accuracy by;
 - aligning (and tracking) computers real-time-clock and datasample clock to prevent this jitter reducing time-stamp accuracy

Hands-on 3: Sequenced Sentences

Experiment Task

- display set of sentences on the screen where every second 1 more character gets added to the sentence
- pause for 5 seconds between sentences (and/or wait for key press)
- send events for everything that happens

Assignment:


- Make flowchart
- Write code -> test -> debug -> until it works :-)
- Start from runSentences_skel.m

Useful Functions:


- sleepSec(time)
 - cause matlab to sleep for the indicated duration in seconds (more accurate than 'pause').



Flowchart: sequenced sentences



Events and processing functions





Hands-on 4: ERP Viewer

Experiment Task

- In 5 sequences of 10 seconds:
 - Every 1 seconds: either randomly display or don't display a cross (+) on the screen for 200ms
- Display a 'Press key to continue string' between sequences, and wait for key press to move to the next sequence
- For every 'stimulus event', i.e. point when the '+' could have been displayed, record 600ms of data annotated with whether it was a '+' or not
- Every time you get some data, compute an average of the EEG data for that type of stimulus, i.e. + or no-+, and display the resulting averages as a multi-plot on the screen
- N.B. You will need a separate signal processing process to: get the data, compute the ERP and display the results!

Assignment:

- Make flowchart for each of the processes, i.e. stimulus, and signalProcessing
- For the expt-control & simulus presentation start from : runStimulus-skel.m
- For the signalProcessing & results genration use : runSigProc-skel.m



Useful Functions:

[data,devents,state]=buffer_waitData(host,port,state,...
'startSet',startEvts,'trlen_samp',samp,'exitSet',exitEvts)

- for all events matching startEvts record samp samples of data
- until an event matching exitEvts is generated
- startEvts and exitEvts specify the events to match in the format:
 - type event type has this value
 - {'type' val} event has type==type and value==val
 - {{'type1' "type2'}} event has type == 'type1' or 'type2'
 - {{'type1' 'type2'} {val1 val2}} event has type == 'type1' or 'type2' and value== val1 or val12
- return the matched event structure(s) in devents and corrospending data in data
 - Data is a vector of structures. data.buf = [nChannels x nSamples] raw EEG data
- state is the match state, used to identify which events have been processed between function calls
- N.B. ExitEvts has the special event type 'data' which returns as soon as the data is available for the first matched startEvt



Summary

- BCI can be broken into 4 processes: dataacquisation, experimental control, signal processing, and stimulus presentation
- buffer_bci framework : uses buffer events as a blackboard for inter-process communication